

HeatProbe Unit and Heat Probe


HeatProbe[™]

More heating power and selectable coagulation levels



OLYMPUS HPU-20 Ensures Reliable Hemostasis

HeatProbe™



REPROCESSING

**Improved reprocessing capability
minimizes
cross-contamination risk**



OPERABILITY

**Two-pedal footswitch for easy
operation**



COAGULATION

**Faster and more reliable — this unique heat-generating device
delivers optimum temperatures
and allows you to choose from six coagulation levels**

With the CD-110U/120U HeatProbe, hemostasis is achieved by a special heat-generating device built into the tip of the probe. In combination with the HPU-20 HeatProbe Unit, which provides higher coagulation heat generating power than those of conventional models, this enables a sure coagulation. You can also fine-tune the system according to the location and

amount of bleeding by selecting one of six different heat energy levels, ranging from 5J to 30J. By adjusting the amount of heat energy according to the level of hemorrhaging, you'll be able to perform more accurate hemostatic treatment. Jet irrigation from the nozzle at the tip washes the site to facilitate positioning, while the tip's special coating prevents sticking to the mucosa.

EFFICIENCY

**Leading-edge thermal control
system**

To reduce the risk of electric shock and heat injury, high-frequency current is not used by the heat probe. Moreover, because the new advanced thermal control system achieves rapid heating and cooling, critical hemostatic treatment can be performed almost instantly.

COMPACT, LIGHTWEIGHT & COST EFFECTIVE

Portable and cost-effective

This lightweight unit is barely two-thirds the size of a conventional model, making it highly portable and especially useful for emergency hemostasis. Above all, this unit is attractively priced to offer a level of performance that's hard to beat.



HeatProbe Unit / Heat Probe

Specifications

■HeatProbe™ Unit (HPU-20)

| | | | |
|-------------------------------|-----------------------|--|----------------------|
| Power supply | Voltage | 120 V AC | 220, 230 or 240 V AC |
| | Input current | 1 A | 0.5 A |
| | Voltage fluctuation | Within $\pm 10\%$ | |
| | Frequency | 50/60 Hz | |
| | Frequency fluctuation | Within ± 1 Hz | |
| Size | Dimensions | 150 (W) \times 170 (H) \times 300 (D) mm | |
| | Weight | 5.9 kg | |
| Compatible Heat Probe™ | | CD-110U:Scope Channel | 3.7 mm |
| | | CD-120U:Scope Channel | 2.8 mm |
| Coagulation | Energy setting | Can be set in 6 steps (5/10/15/20/25/30Joules) using the Energy Setting Switches on the front panel. | |
| | Output | Heat is generated by pressing the “COAG” pedal of the Foot Switch (MAJ-528). Once the Foot Switch pedal is pressed, heat is generated until the set energy value has been supplied. | |
| Water feeding | Water feed setting | Can be set in 4 steps (“High” “.” “Low” “OFF”) using the Water Feed Setting Switches on the front panel. | |
| | Output | Water is fed by pressing the “WASH” pedal of the Foot Switch (MAJ-528). | |

■HeatProbe™

| | | |
|---|-------------|-------------|
| Model | CD-110U | CD-120U |
| Max. insertion portion diameter | 3.4 mm | 2.7 mm |
| Working length | 2300 mm | |
| Channel diameter | 3.7 mm min. | 2.8 mm min. |
| Compatible OLYMPUS HeatProbe™ unit | HPU-20 | |
| Compatible OLYMPUS irrigation pump | MAJ-527 | |

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.



OLYMPUS MEDICAL SYSTEMS CORP.
Shinjuku Monolith, 2-3-1 Nishi-Shinjuku, Shinjuku-ku, Tokyo 163-0914, Japan

For a complete listing of sales and distribution locations visit:
www.olympus.com