

Your Vision, Our Future



Two-in-One

Video System & Light Source Integration



VIDEO SYSTEM

Space-Saving Design

The integrated light source and video processor make the OTV-SI the imaging platform of choice for clinics of all sizes.



Contour Enhancement

Contour enhancement defines contours in the


Contour Enhanced Image

image for clear observation.

Unprecedented Image Quality

To meet the needs for higher quality imaging in ENT

office examinations, we offer the OTV-SI Office Video Endoscopy System that delivers the same true color, high-resolution video imaging.



ENF-VQ Image

Multipurpose Interface for Fiber/Rigid Endoscopes

By connecting a camera head*, you can use your existing fiber/rigid endoscopes.



Normal Image

Exceptional Imaging

Only a Videoscope Can Offer

RHINO-LARYNGO VIDEOSCOPE

ENF-VQ (for diagnosis)

High-Resolution Images

Featuring an all-new CCD that delivers high-resolution images, the ENF-VQ may help detect minute details that would not have been previously detected.

Full-Screen Images : Clear and Bright In Every Corner

Full-screen images make observation much easier. Thanks to a new tip design, the ENF-VQ provides full-screen images that are clear and bright even in the corners.

RHINO-LARYNGO VIDEOSCOPE ENF-V2 (for diagnosis)

Ultra-Slim 3.2 mm Diameter

The ENF-V2 not only delivers some of the brightest, clearest images available from a videoscope, it boasts a revolutionary 3.2 mm ultra-slim diameter. Even thinner than conventional routine fiberscopes, this slim design helps to facilitate insertion while enabling clear observation.

RHINO-LARYNGO VIDEOSCOPE




One Scope Enables Both Observation with Clearer Views and Effective Treatment

- Suction mechanism lets you keep a clear view with easy push-button switching.
- Twenty different types of accessories are available, giving you a wide range of instruments to choose from.
- Compatible with both laser/electrosurgery treatments.


2.0 mm Close Focus Capability


Close focus allows the endoscope tip to come as close as 2.0 mm without electronic magnification, enabling you to observe small details in the mucosa.

All Olympus ENT videoscopes offer outstanding image quality and digital enhancement.










System Chart



Specifications -----

ENF-VQ

Video Sys	stem OTV-SI		
Power Supply	Voltage	100 - 120 / 220 - 240 V AC	Illumination
	Voltage Fluctuation	Within 10%	
	Frequency	50/60 Hz	
	Frequency Fluctuation	Within ±1 Hz	
	Electric Current Consumption	1.2 A	
	Fuse Rating	4 A, 250 V	
	Fuse Size	6.35 x 31.8 mm	
Size	Dimensions	295(W) x 171(H) x 434(D) mm	
	Weight	12 kg	

n	Switching of Lamp	Manual	
	Lamp	150 W halogen lamp (MD-151)	
Lamp Life		Approx. 40 hours of continuous use.	
		When used intermittently, the lamp	
		life may vary slightly.	
	Ignition	Switching Regulator	
	Light Output	Light-path diaphragm control	
	Adjustment		
	Cooling	Forced-air cooling	


Flexible Videoscope

		ENF TYPE VQ	ENF TYPE V2	ENF TYPE VT2
Optical System	Field of View	90°	90°	90°
	Depth of Field	5 - 50 mm	5 - 50 mm	2 - 40 mm
Insertion Tube	Distal End Outer Diameter	3.9 mm	3.2 mm	4.8 mm
	Flexible Tube Outer Diameter	3.6 mm	3.4 mm	4.9 mm
	Working Length	300 mm	300 mm	365 mm
Bending Section	Angulation Range	Up/Down 130°	Up/Down 130°	Up/Down 130°
Total Length		528 mm	528 mm	635 mm
Channel Inner Diameter		-	-	2.0 mm

Note: Product lineup may differ in some countries.



ENF-VT2



ENF-V2

OTV-S7H-N OTV-S7H-1N

Camera Head



Specifications

Camera Head Dimensions 27 x 37 mm (from mount surface) Weight 40 g (excluding cable) Cable 5.2 mm x 4 m Cleaning/Disinfection Immersible in disinfectant solution with a video adapter attached. Sterilization Ethylene oxide gas sterilization possible

•Number of remote control switches: 2

Video Adapter



Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.



OLYMPUS MEDICAL SYSTEMS CORP.

Shinjuku Monolith, 2-3-1 Nishi-Shinjuku, Shinjuku-ku, Tokyo 163-0914, Japan

For a complete listing of sales and distribution locations visit: www.olympus.com