

Single Use Biopsy Valve

More Watertight, More Maneuverable



SINGLE USE BIOPSY VALVE (Sterile)

Single-use design for more safety and convenience

Designed by world-leading endoscope manufacturer OLYMPUS, the MAJ-1555 Single Use Biopsy Valve, fully compatible with your OLYMPUS endoscopes, reduces the risk of cross-contamination and ensures optimal device manipulation. Thanks to an improved new design, the MAJ-1555 prevents body fluid from flowing out of the scope port when a device is inserted during any endoscopic procedures. (ex. Bile juice, blood, and water)



Smoother device manipulation

In spite of enhanced watertightness, the outer slit is longer, preventing the valve from splitting at either end of the slit during device insertion. In addition, the material composition of the new cap decreases the likelihood of the slit closing up during use.



Compatible with a wide range of EndoTherapy devices

Innovative, watertight structure

A new, dual-chamber structure has a smaller inner hole located an appropriate distance from the outer slit. This ensures a better fit for inserted devices, including syringes.



Dual chamber prevents leakage of body fluid out of the scope port. (ex. Bile juice, blood and water)

Secure mounting on the biopsy port

With housing designed to fit securely on the scope's biopsy port, the valve is less likely to loosen or fall off during a procedure. The firm fit helps prevent leakage during syringe injection.



Secure fit prevents valve from loosening or falling off

Easy-to-handle package

Each MAJ-1555 comes in a large, easy-open, sterile package to prevent cross-contamination. Thanks to its larger size, the package can be opened effortlessly during a procedure without fear of the product falling out of its sterile sleeve.



Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

OLYMPUS[®]

OLYMPUS MEDICAL SYSTEMS CORP.
Shinjuku Monolith, 2-3-1 Nishi-Shinjuku, Shinjuku-ku, Tokyo 163-0914, Japan

For a complete listing of sales and distribution locations visit:
www.olympus.com