


WAVICIDE-01[®]

- Ready To Use
- No Activation Required
- 30-Day Re-Use
- Room Temperature
- Safe With Endoscopes



MCC

MEDICAL CHEMICAL CORPORATION

SINCE 1954

WAVICIDE-01®

STERILIZATION/ DISINFECTION SOLUTION	READY TO USE- NO MIXING	STERILIZATION TIME AND TEMPERATURE	NO. OF DAYS FOR STERILIZATION RE-USE	NO. OF DAYS FOR HIGH LEVEL DISINFECTION RE-USE	HIGH-LEVEL DISINFECTION TIME AND TEM- PERATURE	AFTER ACTIVA- TION, SHELF LIFE OF UNUSED PRODUCT
WAVICIDE-01® 2.65% glutaraldehyde	YES	10 Hrs. 22°C	30 Days 22°C	30 Days 22°C	45 Min. 22°C	1 Year
CIDEX® 2.5% glutaraldehyde	NO	10 Hrs. @ 25°C	14 days @ 25°C	14 days @ 25°C	45 min. @ 25°C	14 Days
CIDEX 7® 2.5% glutaraldehyde	NO	10 Hrs. @ 25°C	28 days @ 25°C	28 days @ 25°C	90 min. @ 25°C	28 Days
MetriCide® 2.6% glutaraldehyde	NO	10 Hrs. @ 25°C	14 days @ 25°C	14 days @ 25°C	45 min. @ 25°C	14 Days
MetriCide 28® 2.6% glutaraldehyde	NO	10 Hrs. @ 25°C	28 days @ 25°C	28 days @ 25°C	90 min. @ 25°C	28 Days
CIDEX OPA® 0.55% ortho-Phthalaldehyde	YES	NOT A STERILANT	28 days @ 25°C	14 days @ 20°C	12 min. @ 20°C	28 Days

WAVICIDE-01® AND RELATED PRODUCTS



<u>Description</u>	<u>Catalog #</u>	<u>Packaging</u>
Wavicide-01®	0104	1-gal.
Wavicide-01®	0112	1-qt.
Mec Indicator Strips	0306	50 strips/bottle
Wavizyme® (Concentrate Enzymatic Cleaner)	3304 3312 3336	1-gal. 1-qt. 1-oz. (36/case)
Medizime® (Low-Foam Enzymatic Cleaner)	3304LF 3312LF 3336LF	1-gal. 1-qt. 1-oz. (36/case)



MEDICAL CHEMICAL CORPORATION

19430 Van Ness Avenue • Torrance, CA 90501 • Tel (800) 424-9394 • Fax (310) 787-4464

VISIT OUR WEBSITE AT www.med-chem.com