

# ENDOSCOPE REPROCESSOR OER-Smart

## OER-Smart Related Accessories



### Acecide 1,065ml

High-level Disinfectant

Peracetic acid based liquid chemical germicide designated for OER-Smart.



### Acecide Test Strips

Test strip for Acecide



### EndoRapid 980ml

Alkaline Detergent

Formulated specifically for use with Olympus endoscope reprocessors.



### Gas Filter: MAJ-822

Charcoal filter absorbs disinfectant vapors.

### Air Filter: MAJ-823

Eliminates contaminants suspended in air.

### Water Filter: MAJ-824\*

Removes microorganisms larger than 0.2µm.

\* Water filter, MAJ-2318, may be available in some areas.



### Connecting Tubes

For feeding water or solution into the channels of scopes.



### Connector Hanger: MAJ-865

The connector hanger holds the connector portion of the endoscope and can be mounted on either side of the OER-Smart.

# OLYMPUS

ENDOSCOPE REPROCESSOR

## OER-Smart

Evolved Technology, Proven Reliability



Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

## OLYMPUS

OLYMPUS MEDICAL SYSTEMS CORP.

Shinjuku Monolith, 2-3-1 Nishi-Shinjuku, Shinjuku-ku, Tokyo 163-0914, Japan  
www.olympus-global.com

Global Professional Education website  
www.olympusprofed.com



N8601895-082021