

OLYMPUS[®]


Your Vision, Our Future

Endo**Therapy**[™]

Single Use Hot Biopsy Forceps

EndoJaw Hot[™]

The power of EndoJaw[™] now available for hot biopsy



Single Use Hot Biopsy Forceps with Smooth Insertion Capability and Easy Tissue Targeting

Product Features

Smooth insertion capability

The slim outer diameter combined with the proprietary sheath design and improved flexibility ensures smooth insertion through the endoscope channel.

Specially designed to minimize damage to the scope

As an endoscope manufacturer, OLYMPUS has studied how to minimize the damage to the endoscopes without sacrificing the forceps' performance. Our high quality biopsy forceps will not damage our endoscope biopsy channels.

Stainless steel construction

The biopsy cups are made of stainless steel, which is a material of choice for precision medical devices.

Easy tissue targeting

The forceps cups pivot to allow easier access to the tissue site from an oblique angle.

Full product line-up

Choose either the traditional smooth, oval cup forceps or alligator jaw forceps without needle. Whichever choice you make, all OLYMPUS forceps are made with the quality you can trust and an unconditional guarantee of your satisfaction.

Specifications

Model	Feature	Quantity	Working length	Sheath color	Minimum channel size
FD-210U	Alligator jaw-step	5/box	2300 mm	Blue	2.8 mm
FD-230U	Standard oval	5/box	2300 mm	Blue	2.8 mm



Alligator jaw-step



Standard oval

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.